


DDDDK Tag Antibody / FLAG Tag [clone RM399] (R20415)

Catalog No.	Formulation	Size
R20415-100UG	1 mg/ml in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	DDDDK
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM399
Purity	Protein A purified from animal origin-free supernatant
Applications	Western Blot : 0.001ug/ml-0.1ug/ml
Limitations	This recombinant DDDDK Tag antibody is available for research use only.



Western blot testing of 293T cells (1) transfected or (2) non-transfected with a DNA construct encoding DDDDK tag fusion protein, using recombinant DDDDK Tag antibody at 0.001ug/ml.

Description

This antibody reacts to recombinant proteins containing a DDDDK fused to either the amino or carboxy terminus. No cross reactivity with any endogenous protein.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant DDDDK Tag antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Mixture of a peptide with DDDDK-Tag at the N-terminal and another peptide with DDDDK-Tag at the C-terminal was used

as the immunogen for the DDDDK Tag antibody.

Storage

Store the recombinant DDDDK Tag antibody at -20oC.